



**Roll-to-roll Lab**  
**Syringe Pump**  
**Glove-box**

- w Substrate 110 cm long, up to 20 cm wide.
- w Substrate speed from 5 to 300 cm/min.
- w Electronically controlled substrate tension.
- w Wi-Fi communication allow control by Tablet (Android and iOS) or external computer.
- w Camera analysis module.



For further information, contact us.

[www.developnow.com.br/weler](http://www.developnow.com.br/weler)  
[contato@developnow.com.br](mailto:contato@developnow.com.br)  
São Carlos - Brazil



- w Modular equipment.
- w Easy tool exchange.
- w 6 deposition techniques and 2 drying units.


### Modular Coating Tools




Slot-die



Knife




Wire-bar




Spray

### Modular Printing Tools



Flexography



Rotogravure


### Modular Drying Unit



Oven



Hot plate



For further information, contact us.


[www.developnow.com.br/weler](http://www.developnow.com.br/weler)  
[contato@developnow.com.br](mailto:contato@developnow.com.br)  
São Carlos - Brazil

## Roll-to-roll Lab Syringe Pump Glove-box

- w Compact and lightweight:
  - 85 x 50 x 45 cm (L x W x H)
  - About 40 kg.
- w Drying temperature controlled by PID:
  - Up to 200 °C.
- w Integrated syringe pump.



- w Tablet control of parameters.



For further information, contact us.

[www.developnow.com.br/weler](http://www.developnow.com.br/weler)  
[contato@developnow.com.br](mailto:contato@developnow.com.br)  
São Carlos - Brazil





## Roll-to-roll Lab Syringe Pump Glove-box

- w Compact syringe pump for confined space applications:
  - 195 x 72 x 87 mm (L x W x H)
- w All metal construction, ensuring better strength and finish.
- w Flows between 2.08  $\mu\text{L/hr}$  and 33,3 mL/min (1999 mL/hr) depending on the syringe.



- w One syringe, from smallest size up to 60 mL.
- w Independent control system allows the operation of the pump at distance.
- w Reverse-direction latch allows infuses and withdraws.

\* Not For Clinical Use On Humans

For further information, contact us.

[www.developnow.com.br/weler](http://www.developnow.com.br/weler)  
[contato@developnow.com.br](mailto:contato@developnow.com.br)  
São Carlos - Brazil




## Roll-to-roll Lab Syringe Pump Glove-box

- w Glove-box in stainless steel guarantees better mechanical and solvents resistance.
- w Fully removable side to insert large equipment inside the main chamber.
- w Plugs and electrical cables inside:
  - 100/240 V plugs.
  - Connector and cable for data.
- w LED lighting.
- w Two valves, remove the air and insufflate the working gas.



- w Dimensions: 105 x 60 x 60 cm. (special dimensions on request)
- w High quality Neoprene gloves.
- w Vacuum pump and nitrogen cylinder can be supplied. (optional)
- w Antechamber with 3-way control (nitrogen / closed / vacuum) and individual gauge.



For further information, contact us.

[www.developnow.com.br/weler](http://www.developnow.com.br/weler)  
[contato@developnow.com.br](mailto:contato@developnow.com.br)  
São Carlos - Brazil